

Available on FP10/GP10 - PLEASE COMPLETE ALL SECTIONS
TO BE COMPLETED BY THE HEALTHCARE PROFESSIONAL AND GIVEN TO YOUR PRESCRIBER

1 Patient Name
Doctor/Clinician Clinic
.....
Pharmacist
Date

2 AVAILABLE STYLES:

	BELOW KNEE		THIGH LENGTH	
	Open Toe	Closed Toe	Open Toe	Closed Toe
Class I	Beige	Beige or Black	Beige	Beige or Black
Class II	Beige	Beige or Black	Beige	Beige or Black
Class III	Beige	Beige	Beige	Beige

3 COMPRESSION CLASS:

1 (14-17 mmHg) 2 (18-24 mmHg) 3 (25-35 mmHg)

STYLE:

Below Knee Thigh Length

FOOT:

Closed Toe Open Toe

4 QUANTITY (PAIRS): Beige Black
(Class I & II Closed Toe only)

5 MEASUREMENTS:

MEASUREMENT TIPS:


- Measurements should be taken with the person standing, feet flat on the floor. If your patient is unable to stand, you can take the measurements with the person seated, and feet flat on the floor.
- For foot measurement: stand on a plain sheet of paper, mark at the toe and at the heel; remove the foot and measure between the two points; this will provide a more accurate reading.
- When possible, measurements should be taken first thing in the morning. ⁽¹⁾
- For thigh length, take the measurement while the person is standing.

1. NICE Clinical Knowledge Summary - Compression Stockings (2008, last update 2022)

All measurements are in centimetres

Mid thigh (A) measurement is only required for thigh length stockings

Foot (D) measurement is only required for closed toe stockings



	A: Mid Thigh	B: Calf	C: Ankle	D: Foot	Size Required
Small	41 - 54	32 - 38	19 - 25.5	20 - 26	
Medium	44 - 58	34 - 40.5	21 - 27.5	22 - 28	
Large	47 - 61	36 - 43.5	23 - 29.5	24 - 31	
X Large	50 - 64	38 - 47	25 - 32	26 - 33	

COMMENTS:

NEED HELP? Contact Urgo Medical for CPD accredited training. **T: 0330 128 0898**